



SIGMAX

GENERAL CATALOGUE



Contact Information:

Company Name:

Email:

Phone:

Fax:

NIPPON SIGMAX CO., LTD.

33rd Floor Shinjuku Oak Tower 6-8-1

Nishi-shinjuku, Shinjuku-ku, Tokyo 163-6033 JAPAN

www.sigmax-med.jp/en

BRACES & SUPPORTS

1

BRACES & SUPPORTS

- ▶ Ankle P3
- ▶ Knee P10
- ▶ Back P17
- ▶ Rib P21
- ▶ Neck P23
- ▶ Others P25

2

BANDAGE PRODUCTS

- ▶ Adhesive Elastic Bandage P33
- ▶ Elastic Bandage P34
- ▶ Net & Tubular Bandage P35

3

CASTING PRODUCTS

- ▶ Aluminum Splints P37
- ▶ Soft Splints P37
- ▶ Others P38

▶ Back



P17_20

▶ Neck



P23_24

▶ Rib



P21_22

▶ Knee



P10_16

▶ Ankle




P3_9

▶ Others



P25_31





Designed to fit well,
be comfortable and provide
excellent stability.

ANKLE

FACILIAID ANKLE

[Indication] ● Ankle joint support



Comfortable Material

Soft, elastic material is used for this brace. It is made of an antimicrobial finished fabric and washable, so it can be reused hygienically.



Moderate Compression

The wrap uses figure 8 taping techniques to stabilize and protect the ankle joint. The elastic material and Velcro fastener make the wrap compatible with a wide range of ankle joint sizes. The amount of compression and fit can be adjusted.



Two Colors Available

Beige and black colors are available.



Easy Application

The end of the wrap is loop-shaped, allowing this to be easily wrapped around the ankle.



How to Apply

STEP 1

Pass the foot through



STEP 2

Wrap



STEP 3

Fasten



For both the right and left feet

Item Number	Size	Color	Shoe Size	Qty per Box
302501	One Size	Beige	20.0 - 27.5 cm	5 Bags
302502	One Size	Black	20.0 - 27.5 cm	5 Bags

FACILIAID ANKLE SOFT

- [Indication]**
- Medial, lateral ankle instability
 - Ankle joint protection

FACILIAID ANKLE SOFT is an ankle brace that provides excellent mobility to support the user's lifestyle and comfortable fit while providing support for the ankle joint.

It is easy to put on and take off, and can be used in recovery and to prevent further occurrence of the injury.



Breathable Material

Reduces thickness by using thin, elastic STRAFLEX™ as the core material. It can be worn with business shoes such as pumps. The breathable material prevents clamminess when worn for long periods.



*STRAFLEX™ has excellent elasticity.

Figure 8 Strap

Supports the ankle with a structure based on the figure 8 taping technique. The hook-and-loop fastener allows to attach any part of the soft fabric, allowing fit and the amount of compression to be adjusted freely.



Open-Back Design

The open-back design makes the brace easy to put on. As the foot is inserted from the back, you can simply place your foot in the accurate position.



Antimicrobial and Odor-Resistant

Adds antimicrobial and odor-resistant finishing to the fabric to improve hygiene and make the brace more comfortable.

For both right and left feet

Item Number	Size	Shoe Size
F303502	M	19 - 25 cm
F303503	L	25 - 31 cm

*If the shoe size falls between two sizes, choose the larger one.

EXAID ANKLE 3

- [Indication]**
- Lateral ankle instability
 - Ankle sprain

EXAID ANKLE 3 provides lateral ankle support.

The original stay, which is designed to fit the ankle joint, prevents inversion and improves comfort.

The brace is also made from a hygienic and comfortable fabric.



Original Medial Stay

The built-in original stay is designed to prevent inversion without obstructing natural daily life movement. The stay is covered to reduce discomfort when wearing shoes.




Elastic Strap

The elastic strap enhances support provided by the original stay by wrapping around the plantar and dorsal of the foot from the inside.



Breathable Material

Reduces the thickness by using thin, elastic STRAFLEX™ as the core material. It can be worn with business shoes such as pumps. The breathable material prevents clamminess when worn for long periods.



*STRAFLEX™ has excellent elasticity.

Antimicrobial and Odor-Resistant

Adds antimicrobial and odor-resistant finishing to the fabric to improve hygiene and make the brace more comfortable.

Open-Back Design

The open-back design makes the brace easy to put on. As the foot is inserted from the back, you can simply place your foot in the accurate position.



Item Number	Size	Shoe Size
F452301	Right S	19 - 22 cm
F452302	Right M	22 - 25 cm
F452303	Right L	25 - 28 cm
F452304	Right XL	28 - 31 cm

*If the shoe size falls between two sizes, choose the larger one.

Item Number	Size	Shoe Size
F452311	Left S	19 - 22 cm
F452312	Left M	22 - 25 cm
F452313	Left L	25 - 28 cm
F452314	Left XL	28 - 31 cm

*Wear the product over socks according to the instructions.

ANKLE BRACES & SUPPORTS
KNEE BRACES & SUPPORTS
BACK BRACES & SUPPORTS
RIB BRACES & SUPPORTS
NECK BRACES & SUPPORTS
OTHERS
NURSING CARE PRODUCTS
ADHESIVE ELASTIC BANDAGE
ELASTIC BANDAGE
NET BANDAGE/TUBULAR BANDAGE
CASTING PRODUCTS
ALUMINIUM SPLINTS
SOFT SPLINTS
OTHERS

ANKLE BRACES & SUPPORTS
KNEE BRACES & SUPPORTS
BACK BRACES & SUPPORTS
RIB BRACES & SUPPORTS
NECK BRACES & SUPPORTS
OTHERS
NURSING CARE PRODUCTS
ADHESIVE ELASTIC BANDAGE
ELASTIC BANDAGE
NET BANDAGE/TUBULAR BANDAGE
CASTING PRODUCTS
ALUMINIUM SPLINTS
SOFT SPLINTS
OTHERS

EXAID ANKLE 6

- [Indication]**
- Medial, lateral and anterior ankle instability
 - Ankle sprain

EXAID ANKLE 6 provides medial and lateral ankle support. You can adjust the stability according to the situation. It has many innovative design features that improve mobility and comfort.



Original Rigid Guard


The original medial and lateral guards have been designed separately, focusing on the difference in position of the medial and lateral malleolus to enhance stabilization and fit. The joint plate is built in at the ankle part to make it possible to walk naturally.



The well-fitting 3D structure and the lightweight guard improve mobility and make it easier to wear shoes over the brace.



The connection part of the medial/lateral guards is flexible to reduce discomfort.




Open-Back Design

The open-back design makes the brace easy to apply. As the foot is inserted from the back, you can simply place your foot in the accurate position.



Figure 8 Strap


The strap crosses at the front of the ankle bone to form a figure of 8, providing anterior, medial and lateral support based on the figure 8 taping techniques.



Item Number	Size	Shoe Size
F452601	Right S	19 - 22 cm
F452602	Right M	22 - 25 cm
F452603	Right L	25 - 28 cm
F452604	Right XL	28 - 31 cm

* If the shoe size falls between two sizes, choose the larger one.

Design of the Original Guard



- Curved line that becomes wider at the top reduces discomfort.
- V-shaped reinforcing rib enhances support and strength.
- The connection part of the guard has been widened to be compatible with various calf sizes.
- The rigid parts such as plates and rivets are flatter to reduce discomfort, and have been made more durable.

Item Number	Size	Shoe Size
F452611	Left S	19 - 22 cm
F452612	Left M	22 - 25 cm
F452613	Left L	25 - 28 cm
F452614	Left XL	28 - 31 cm

*vWear the product over socks according to the instructions.

EXAID ANKLE FO

- [Indication]**
- Medial, lateral and anterior ankle instability
 - Ankle sprain

EXAID ANKLE FO has an open-front design that wraps around ankles of any size and FRTP guard that provides strong support. These features make the brace fit as if it has been custom-made, and provide targeted support.



Open-Front Design

The open-front design makes the brace easy to apply and fit the ankle shape.



FRTP Guard

The FRTP (fiberglass reinforced thermoplastic) guard has been designed based on the average of anatomical data collected by the 3D measurements for each size to reduce discomfort. When necessary, fitting can be further improved by heating the guard and molding it to fit the foot. The FRTP guard is also removable, so can be used or not used according to the ankle condition.





Figure 8 Strap

The strap crosses at the front of the ankle bone to form a figure of 8, providing anterior, medial and lateral support based on the figure 8 taping techniques.




Heel Strap

The heel strap enhances support by connecting to the FRTP guards on both sides.



Breathable Material

The plantar section of the brace is made of thin, elastic STRAFLEX™. The breathable material prevents clamminess when worn for long periods.



Item Number	Size	Shoe Size
F452801	Right S	19 - 22 cm
F452802	Right M	22 - 25 cm
F452803	Right L	25 - 28 cm
F452804	Right XL	28 - 31 cm

* If the shoe size falls between two sizes, choose the larger one.

Item Number	Size	Shoe Size
F452811	Left S	19 - 22 cm
F452812	Left M	22 - 25 cm
F452813	Left L	25 - 28 cm
F452814	Left XL	28 - 31 cm

*Wear the product over socks according to the instructions.

FACILIAID ACHILLES

[Indication] ● Protection of the Achilles tendon




3D-Cut Fabric

The 3D-cut body and outer strap provide a secure fit.



Achilles Pad

The Achilles pad supports the tendon gently and provides stabilization.



Non-Elastic Straps

The outer strap is made from non-elastic material to restrict ankle dorsiflexion that stretches the Achilles tendon.



Heel Lock

The taping (Heel Lock) feature on the outer strap increases heel stabilization. Stabilizing the calcaneus reduces stress on the Achilles tendon.



*In order to clearly demonstrate how the product fits, the pictures show the brace being worn on bare foot. It should be worn over socks during actual use.

For both right and left feet

Item Number	Size	Shoe Size
F302202	M	21 - 25 cm
F302203	L	25 - 29 cm

*If the shoe size falls between two sizes, choose the larger one.

Wearing the brace on the right leg



■ Front ■ Medial ■ Back ■ Lateral



K N E E



**Our ergonomic designs
are user-friendly and ideal
for use over long periods of time.**

K N E E

EXAID KNEE MCL & LCL

[Indication]

- Provides medial support for MCL (Medial Collateral Ligament) injuries
- Provides lateral support for LCL (Lateral Collateral Ligament) injuries

(Features)

- Antimicrobial, Odor-Resistant
- 4-Way Stretch Mesh Material
- Binder Tape
- Aluminum Offset Joints



Aluminum Joints

The aluminum offset joints improve lateral /medial stability.



3-Point Support

The multi-function strap provides 3-point support. This reduces medial (for MCL injuries) and lateral (for LCL injuries) instability.



Fastening Strap

Non-elastic tape inside the strap provides firm, constant support. It reduces anterior instability in the tibia and is suitable for combined MCL and LCL injuries.



Easy Application

The brace can be put on, taken off and adjusted easily at the front side without any need to reach around to the back of the knee.



*Reverse the right and left sides of the product for LCL injuries.

Item Number	Size	Thigh Circumference
F363001	S	34 - 39 cm
F363002	M	39 - 44 cm
F363003	L	44 - 50 cm
F363004	XL	50 - 56 cm
F363011	S	34 - 39 cm
F363012	M	39 - 44 cm
F363013	L	44 - 50 cm
F363014	XL	50 - 56 cm

[Taking measurements and choosing the correct size] Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



Item Number	Size	Thigh Circumference
F363101	S	34 - 39 cm
F363102	M	39 - 44 cm
F363103	L	44 - 50 cm
F363104	XL	50 - 56 cm

[Taking measurements and choosing the correct size] Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.

EXAID KNEE ACL

[Indication]

- Provides anterior tibia support for ACL (Anterior Cruciate Ligament) injuries

(Features)

- Antimicrobial, Odor-Resistant
- 4-Way Stretch Mesh Material
- Binder Tape
- Aluminum Offset Joints

Aluminum Joints




The aluminum offset joints improve lateral /medial stability.

4-Point Support



The multi-function strap, fastening strap and elastic strap provide 4-point support. This reduces anterior instability in the tibia.

Fastening Strap



Non-elastic tape is built in to the fastening strap provides firm, constant support.

Easy Application

The brace can be put on, taken off and adjusted easily at the front side without any need to reach around to the back of the knee.

EXAID KNEE PCL

[Indication]

- Provides posterior tibia support for PCL (Posterior Cruciate Ligament) injuries

(Features)

- Antimicrobial, Odor-Resistant
- 4-Way Stretch Mesh Material
- Binder Tape
- Aluminum Offset Joints

Aluminum Joints



The aluminum offset joints improve lateral /medial stability.

4-Point Support



The multi-function strap, fastening strap and elastic strap provide 4-point support. This reduces posterior instability in the tibia.

Fastening Strap



Non-elastic tape is built in to the fastening strap provides firm, constant support.

Easy Application

The brace can be put on, taken off and adjusted easily at the front side without any need to reach around to the back of the knee.



Item Number	Size	Thigh Circumference
F363201	S	34 - 39 cm
F363202	M	39 - 44 cm
F363203	L	44 - 50 cm
F363204	XL	50 - 56 cm

[Taking measurements and choosing the correct size] Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



EXAID KNEE Light Sports 2

[Indication] ● Knee support

(Features)

- Antimicrobial, Odor-Resistant
- 4-Way Stretch Mesh Material
- Binder Tape
- Active Resin Stays

Active Resin Stays **Secure Fit** **Open-Front Design**



Provides a secure fit without curling, twisting or bending.

The design of the brace prioritizes stretch, fit and knee mobility.

The open-front design makes the brace easy to apply and remove.

Item Number	Size	Thigh Circumference
F362001	S	34 - 39 cm
F362002	M	39 - 44 cm
F362003	L	44 - 50 cm
F362004	XL	50 - 56 cm

[Taking measurements and choosing the correct size]
Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



EXAID KNEE Patella

[Indication] ● Patella Stabilization
● Patella Subluxation and chondromalacia patella

(Features)

- Antimicrobial, Odor-Resistant
- 4-Way Stretch Thin Material
- Binder Tape
- Active Resin Stays

Active Resin Stays **Patella Pad** **Unique Design**



Provides a secure fit without curling, twisting or bending.

The built-in patella pad reduces lateral deviation of the kneecap. The brace applies traction to the patella pad with the center strap, and stabilizes the pad with the upper and lower straps, allowing the pad to securely fit around the kneecap.

The patella pad does not contact with the skin, which reduces friction.

*The images are for the right knee. If you apply on the left knee, turn the product upside down.

Item Number	Size	Thigh Circumference
F364001	S	34 - 39 cm
F364002	M	39 - 44 cm
F364003	L	44 - 50 cm
F364004	XL	50 - 56 cm

[Taking measurements and choosing the correct size]
Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



EXAID KNEE Light Sports 3

[Indication] ● Knee support
● Medically / Laterally instabilities

(Features)

- Antimicrobial, Odor-Resistant
- 4-way Stretch Mesh Material
- Binder Tape
- Aluminum Offset Joints

Aluminum Joints **Excellent Fit** **Open-Front Design**



The aluminum offset joints improve lateral/ medial stability.

The stretchable material fits the motion of kneecap.

The open-front design makes the brace easy to apply and remove.

Item Number	Size	Thigh Circumference
F362101	S	34 - 39 cm
F362102	M	39 - 44 cm
F362103	L	44 - 50 cm
F362104	XL	50 - 56 cm

[Taking measurements and choosing the correct size]
Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



EXAID KNEE Osgood

[Indication] ● Ideal for reduce Osgood-Schlatter

(Features)

- Antimicrobial, Odor-Resistant
- 4-Way Stretch Thin Material
- Binder Tape
- Active Resin Stays

Osgood Pad **Osgood Band** **Easy Application**



The Osgood pad applies compression to the patellar ligament to relieve pain. The cover for the Osgood pad improves cushioning for a more comfortable fit.

Fasten the Osgood band after fitting the brace, to provide targeted appropriate compression.

The tab makes it easier to hold the brace when apply it.

[Taking measurements and choosing the correct size]
Measure the circumference of the leg 5 cm below the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.

Item Number	Size	Thigh Circumference
F364101	S	29 - 32 cm
F364102	M	32 - 35 cm
F364103	L	35 - 39 cm




EXAID KNEE Short

[Indication] ● Medial / Lateral knee instability

(Features)

- Antimicrobial, Odor-Resistant
- Binder Tape
- Active Resin Stays

Active Resin Stays



Provides a secure fit without curling, twisting or bending.

Mesh Material



Soft, breathable material is used at the back of the knee to prevent clamminess.


Easy Application



The brace has a hole in the front, into which a finger can be hooked to make the brace easier to put on and take off.

Compression Pad

The patella compression pad fits around the kneecap to improve its stability.



Short Length

The brace is shorter to improve mobility and makes it easier to use. It provides knee joint support for those suffering from mild knee osteoarthritis.

Item Number	Size	Thigh Circumference
F361001	S	34 - 38 cm
F361002	M	38 - 42 cm
F361003	L	42 - 46 cm
F361004	XL	46 - 50 cm

[Taking measurements and choosing the correct size]
Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



FACILIAID KNEE

[Indication] ● Knee joint support

(Features)

- 4-Way Stretch Thin Material
- Binder Tape
- Active Resin Stays

Easy Application



The brace has a hole in the front, into which a finger can be hooked to make the brace easier to put on and take off.

Mesh Material



Soft, breathable material is used at the back of the knee to prevent clamminess.

Active Resin Stays



Provides a secure fit without curling, twisting or bending.

Item Number	Size	Thigh Circumference
F302401	S	34 - 38 cm
F302402	M	38 - 44 cm
F302403	L	44 - 50 cm
F302404	XL	50 - 56 cm

[Taking measurements and choosing the correct size]
Measure the circumference of the thigh 10 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.


KNEE SPLINT 〈Extension Type〉 〈Flexion Type〉

[Indication] ● Splinting for knee joint, post knee joint surgery



Adjustable Stays

- Both sides of the knee splint have two tubular pockets (stay covers). The position of the stay can be adjusted according to which pocket it is placed in, the front or the back.
- The stay on the back of the knee is made from strong aluminum. It can be adjusted to suit the shape of the patient's leg.




Comfortable Fit

- Cushioning has been added to the back of the knee splint for improved fit. It also relieves local pressure exerted by the stay at the back.



Cushioning Material



KNEE SPLINT 〈Extension Type〉

Item Number	Size	Length	Thigh Circumference	Thigh Length
301001	S	48.5 cm	33 - 41 cm	26 cm and more
301002	M·L	48.5 cm	41 - 53 cm	26 cm and more
301012	Long M·L	52.0 cm	41 - 53 cm	29 cm and more
301014	Long XL	52.0 cm	53 - 65 cm	29 cm and more
301021	Short S	41.5 cm	33 - 41 cm	19 cm and more
301022	Short M·L	41.5 cm	41 - 53 cm	19 cm and more

KNEE SPLINT 〈Flexion Type〉



Item Number	Size	Length	Thigh Circumference	Thigh Length
301101	S	48.5 cm	33 - 41 cm	26 cm and more
301102	M·L	48.5 cm	41 - 53 cm	26 cm and more
301112	Long M·L	52.0 cm	41 - 53 cm	29 cm and more
301114	Long XL	52.0 cm	53 - 65 cm	26 cm and more


*Select a size based on thigh circumference. Measure the circumference of the thigh 10 cm above the center of the kneecap. If the measurement falls between two sizes, choose the larger one.

*Use the inside leg measurement of the thigh as a reference when choosing a product length. Measure the inside of the leg from the height of the center of the kneecap to the groin.


Breathable Material

- Breathable material at the back of the knee splint to decrease clamminess and make the knee splint more comfortable.





**Complete stability,
mobility, and comfort with
our ergonomic designs.**



BACK

BACK BRACES

MAXBELT me1

Two soft stays are built into the back portion.
Stabilize the back without restricting everyday movement.


Stability ●●●●●
Mobility ●●●●●
Comfort ●●●●●

Light Support

Contains Latex Rubber




Soft Stays Two Soft Stays
Soft stays to enhance support
Soft stays are used to prioritize mobility. They enhance back stability without preventing movement.



Sub-Straps
Sub-straps with pockets
The front of the sub-straps are designed with pockets. A finger can be hooked into this pocket to help to pull the belt, making it ideal for those weak people.



Breathable Material
Two kinds of mesh materials for a comfortable fit
A hard mesh material is used in the back, and an elastic mesh material is used on the sides for improved breathability and increased comfort.



Size	XS	S	M	L	XL	2XL
Item Number	F322100	F322101	F322102	F322103	F322104	F322105
Waist Circumference	55 - 65 cm	65 - 75 cm	75 - 85 cm	85 - 95 cm	95 - 105 cm	105 - 115 cm

MAXBELT me2

Two S-shaped lightweight stays are built into the back portion.
Provides firm support and enhanced back stabilization.
Standard and universal mesh type.


Stability ●●●●●
Mobility ●●●●●
Comfort ●●●●●

Moderate Support

Contains Latex Rubber




S-Shaped Stays Two S-Shaped Stays
S-shaped lightweight stays fit the back
Long S-shaped stays fit the curve of the back and provide strong support. Realized both strong support and weight reduction by changing the thickness of each part of the stay.



Sub-Straps
Sub-straps increase abdominal pressure
The sub-straps enhance compression along the body from the back and increase abdominal pressure to support the back.



Breathable Material
Two kinds of mesh materials for a comfortable fit
A hard mesh material is used in the back, and an elastic mesh material is used on the sides for improved breathability and increased comfort.



Stability
Dual layered structure to enhance stabilization



Size	XS	S	M	L	XL	2XL	3XL
Item Number	F322200	F322201	F322202	F322203	F322204	F322205	F322206
Waist Circumference	55 - 65 cm	65 - 75 cm	75 - 85 cm	85 - 95 cm	95 - 105 cm	105 - 115 cm	115 - 125 cm

*Measure the waist circumference around the navel. If the measurement falls between two sizes, choose the larger one.

MAXBELT me3

3D-cut hard mesh material in the back provides strong support. Two S-shaped lightweight stays for more support and stability.

Stability ●●●●●
 Mobility ●●●●●
 Comfort ●●●●●

Strong Support

Contains Latex Rubber




S-Shaped Stays

Two S-Shaped stays

S-shaped lightweight stays to fit the waist

Long S-shaped stays fit the curve of the back and provide strong support. Realized both strong support and weight reduction by changing the thickness of each part of the stay.



Sub-Straps

Sub-straps increase abdominal pressure


The sub-straps enhance compression along the body from the back and increase abdominal pressure to support the back.



Breathable Material

Two kinds of mesh material for a comfortable fit

A hard mesh material is used in the back, and an elastic mesh material is used on the sides for improved breathability and increased comfort.



3D Structure

3D structure to provide firm back support

A 3D-cut hard mesh material has been used in the back. This provides a good fit and strong support around the back.



Size	XS	S	M	L	XL	2XL	3XL
Item Number	F322300	F322301	F322302	F322303	F322304	F322305	F322306
Waist Circumference	55 - 65 cm	65 - 75 cm	75 - 85 cm	85 - 95 cm	95 - 105cm	105 - 115 cm	115 - 125 cm

MAXBELT me black

Mesh type with sporty black color ideal for daily life. A hard mesh material is used both in the front and the back. Tapered front belt provides comfortable fit.

Stability ●●●●●
 Mobility ●●●●●
 Comfort ●●●●●

Moderate Support




S-Shaped Stays

Two S-Shaped Stays

S-shaped lightweight stays fit the waist

Long S-shaped stays fit the curve of the back and provide strong support. Realized both strong support and weight reduction by changing the thickness of each part of the stay.



Sub-Straps

Sub-straps increase abdominal pressure


The sub-straps enhance compression along the body from the back and increase abdominal pressure to support the back.



Breathable Material

Two kinds of mesh materials provide breath ability and comfort

A hard mesh material is used in the back and front, and an elastic mesh material is used on the sides for improved breathability, increased comfort and a secure fit.



Black Color

Sporty black color ideal for daily life

The sporty black color is useful for everyday use. The non-bulky design allows the brace to be worn under clothes without leaving any noticeable lines.



Size	XS	S	M	L	XL	2XL	3XL
Item Number	F322400	F322401	F322402	F322403	F322404	F322405	F322406
Waist Circumference	55 - 65 cm	65 - 75 cm	75 - 85 cm	85 - 95 cm	95 - 105 cm	105 - 115 cm	115 - 125 cm

*Measure the waist circumference around the navel. If the measurement falls between two sizes, choose the larger one.


MAXBELT E2

Short length type which provides improved comfort. Thin, non-bulky material does not leave noticeable lines when worn under clothes. Mesh material is used to provide breathability and comfort.

Stability ●●●●●
 Mobility ●●●●●
 Comfort ●●●●●

Moderate Support

Contains Latex Rubber



Plastic Stays

Two Plastic Stays

Plastic stays enhance stability

Two flexible plastic stays are built-in to enhance back stability without preventing movement.



Sub-Straps

Sub-straps increase abdominal pressure

The sub-straps enhance compression along the body from the back and increase abdominal pressure to support the back.



Breathable Material

Breathable mesh material

Breathable mesh material provides comfort and prevents clamminess when worn over clothes or underwear.



Easy Application

Easy-to-wear design

Designed to be short for natural movement and decreased discomfort.



Size	XS	S	M	L	XL	2XL
Item Number	F324200	F324201	F324202	F324203	F324204	F324205
Waist Circumference	55 - 65 cm	65 - 75 cm	75 - 85 cm	85 - 95 cm	95 - 105 cm	105 - 115 cm

MAXBELT SHORT

The narrow belt which provides compression to the pelvis and fits tightly without preventing movement.

Stability ●●●●●
 Mobility ●●●●●
 Comfort ●●●●●

Light Support

Narrow Type
 Pelvic Belt



Sub-Straps

Comfortable even over long periods of time


The narrow shape applies compression to the pelvis and is easy-to-wear in daily life.



Comfortable Material

Breathable and comfortable

Breathable material is used to prevent clamminess.



Size	S	M	L	XL	2XL
Item Number	F324101	F324102	F324103	F324104	F324105
Inferior Abdominal Circumference	75 - 83 cm	83 - 91 cm	91 - 99 cm	99 - 107 cm	107 - 115 cm

*Select a size based on inferior abdominal circumference (at the level of the anterior superior iliac spine). Refer to the image on the left for the measurement position.


Well-Fitting & Hygienic – RIB BAND Series



RIB

Breathable fabric is used for the main body, and a soft pile fabric is used inside.

Both ends of the brace have pockets on them, so the brace can be easily pulled tight and secured.



Front of the RIB BAND



RIB BAND



RIB BAND Women



RIB BAND II



RIB BAND me



RIB BAND Fit

The product is designed with breathability and comfort in mind, using an antimicrobial, odor-resistant finishing and mesh material, as it is worn directly on the body.



Back of the RIB BAND

Breathable, elastic fabric is used for the back side. The bottom of this fabric has been narrowed to improve fit.



The shape is designed not to press into the armpits, and provides a secure fit around the upper ribs.

RIB BRACES

[Indication] ● Compression and stabilization of the rib

Standard Type RIB BAND / RIB BAND Women



- Available for various body types with a wide range of sizes. There are bands both for general use and for women.
- Polyurethane foam with an antimicrobial, odor-resistant finishing is used for the main body.



RIB BAND Women is designed not to compress the bust.

Item Number	Size	Chest Circumference*
F313000	XS	55 - 70 cm
F313001	S	70 - 85 cm
F313002	M	85 - 100 cm
F313003	L	100 - 120 cm
F313004	XL	120 - 140 cm
F313005	2XL	140 - 160 cm
F313102	M	60 - 75 cm
F313103	L	75 - 90 cm
F313104	XL	90 - 105cm

*For the general use type, measure the chest circumference. For the women's type, measure the under bust circumference.

Standard Mesh Type RIB BAND me



Item Number	Size	Chest Circumference
F313301	S	70 - 85 cm
F313302	M	85 - 100 cm
F313303	L	100 - 120 cm
F313304	XL	120 - 140 cm

- Breathable and absorbent mesh material is used for the band. The inner part is made of a soft, comfortable texture pile fabric with no stiffness.
- Breathable and elastic mesh material is used in the back side.

Two Straps Regular Type RIB BAND II



Item Number	Size	Chest Circumference
F313201	S	70 - 85 cm
F313202	M	85 - 100 cm
F313203	L	100 - 120 cm
F313204	XL	120 - 140 cm

- The front part is divided into a lower and upper part so compression can be adjusted according to body type and symptoms.
- Polyurethane foam with an antimicrobial, odor-resistant finishing is used for the main body.

Two Straps Mesh Type RIBBAND Fit



Item Number	Size	Chest Circumference
F313401	S	60 - 80 cm
F313402	M	80 - 100 cm
F313403	L	100 - 120 cm
F313404	XL	120 - 140 cm
F313405	2XL	140 - 160 cm

- The front part is divided into a lower and upper parts. So compression can be adjusted according to body type and symptoms.
- Thin mesh fabric is used as the elastic strap to reduce bulkiness.
- V-shaped cut provides more compression at the both upper and lower edges than the center of the support which prevents slipping.

Comfortable & Hygienic

- COLLAR KEEPER Series



As the area around the neck can easily become sweaty, polyurethane foam with an antimicrobial, odor-resistant finish is used for the core.



COLLAR KEEPER



COLLAR KEEPER U

A soft comfortable cotton material is used for the outer cover.

Both the cover and the core are washable. You can always keep them clean.



The removable cover is easy to wash. (Only for COLLAR KEEPER)



Two sheets can be used to adjust the height of the brace, and the curved shape fits the chin and jaw securely.



COLLAR KEEPER me

The double stopper fastening to prevent the brace to come loose.

For safety, the formaldehyde content is lower than permissible levels for baby products.




NECK BRACES

[Indication] ● Stabilization of the cervical region



Standard Type COLLAR KEEPER

- Soft polyurethane foam is used for the inner parts which touch the jaw and neck, and hard polyethylene foam is used for the outer parts to provide both comfort and a secure fit.
- Shape-memory resin is built into the core to help form an easily wearable shape.




Item Number	Size	Neck Circumference	Product Height
F311000	XS	25 - 30 cm	7.5 cm
F311001	S	30 - 35 cm	8.0 cm
F311002	M	35 - 40 cm	8.5 cm
F311003	L	40 - 45 cm	8.0 cm

Contains Latex Rubber



Washable - Cover Type COLLAR KEEPER U

- The brace is designed to fit the shape of the neck. Soft, elastic polyurethane foam is used for the core.
- The brace has good breathability with many air holes.
- Whole brace is washable to keep it clean.




Item Number	Size	Neck Circumference	Product Height
F311101	S	30 - 35 cm	8.0 cm
F311102	M	35 - 40 cm	8.5 cm
F311103	L	40 - 45 cm	8.0 cm




Adjustable Mesh Type COLLAR KEEPER me

- Mesh material provides excellent breathability.
- Two sheets can be moved up and down to adjust the height.
- To provide a secure fit, the top part naturally curves to fit the jaw, and the lower part naturally curves to fit the clavicle.



Item Number	Size	Neck Circumference	Product Height
F311301	S	30 - 35 cm	8.0 - 9.5 cm
F311302	M	35 - 40 cm	8.5 - 10.0 cm
F311303	L	40 - 45 cm	8.0 - 9.5 cm



We have variety of products for each body parts over a wide range of symptoms.



CLAVICLE BRACE

[Indication]

- Stabilization after clavicle fracture repositioning

- Provides moderate compression with thick cushioning material.
- Reduced compression on the armpit.
- The angle of the strap cushion can be adjusted based on shoulder width.
- The double fastening system of an adjuster and hook-and-loop fastener allows for secure fastening.
- Twisting is prevented with a plastic stay in the back pad.
- Wide range of sizes from 2XS to XL.



Item Number	Size	Reference Measurements	
		Men (Chest Circumference)	Women (Under bust)
F15601	2XS	40 - 50 cm	40 - 50 cm
F15602	XS	50 - 60 cm	50 - 60 cm
F15603	S	60 - 70 cm	60 - 65 cm
F15604	M	70 - 85 cm	65 - 80 cm
F15605	L	85 - 115 cm	80 - 115 cm
F15606	XL	115 - 150 cm	115 - 150 cm


*Measure the chest circumference for men, and the underbust circumference for women.

HIP STABILIZER

[Indication]

- Compression of waist and thigh after artificial joint replacement or epidermization

- The two hook-and-loop fasteners on the back and side of the thigh can be used to adjust compression.
- Elastic, well-fitting fabric is used for the main body of the brace.
- The material is comfortable, with excellent absorption and moisture permeability.
- Can be used on either right or left leg.



◀ Two hook-and-loop fasteners make the brace easy to wear and adjust the level of compression.

Item Number	Size	Measurements	
		Hip Circumference	Thigh Circumference
F17001	S	60 - 75 cm	40 - 45 cm
F17002	M	75 - 90 cm	45 - 51 cm
F17003	L	90 - 110 cm	51 - 58 cm
F17004	XL	110 - 130 cm	58 - 65 cm

*Select a size based on hip circumference, and use the thigh circumference as a guide.



WOOL SUPPORTS

[Indication] ● Joint compression and retains heat

- Excellent heat retention with wool fabric.
- Provides compression to joints and keeps them warm.
- Stretches in all directions to gently fit the affected area.
- Made of a blended wool fabric.

Contains Latex Rubber

Item Number	Type	*Measurements	Qty per Bag	Qty per Box
301301	Elbow S	20 - 25 cm	1	6 Bags
301302	Elbow M·Knee S	23 - 30 cm	1	6 Bags
301303	Elbow L·Knee M	30 - 40 cm	1	6 Bags
301304	Knee L	40 - 50 cm	1	6 Bags
301305	Knee M (Long)	30 - 40 cm	1	6 Bags
301306	Palm	—	1	6 Bags
301307	Wrist	17 - 23 cm	2	6 Bags
301308	Ankle	(Shoe Size) 23 - 27 cm	1	6 Bags
301309	Shoulder M	(Chest Circumference) 65 - 75 cm	1	3 Bags
301310	Shoulder L	(Chest Circumference) 75 - 95 cm	1	3 Bags

*For the elbow, measure the arm circumference 10 cm above the center of the elbow. For the knee, measure the leg circumference 12 cm above the center of the kneecap. For the ankle, use shoe size as a guide.

SEAMLESS SUPPORTS

[Indication] ● Joint compression and retains heat

- Economy type with an excellent soft texture.
- Provides compression to joints and keeps them warm.
- Made of a soft material with excellent moisture absorption.
- Available for various uses.

Contains Latex Rubber

Item Number	Type	*Measurements	Qty per Bag	Qty per Box
301401	Wrist	15 - 20 cm	2	6 Bags
301402	Palm	—	1	6 Bags
301403	Elbow	20 - 25 cm	1	6 Bags
301404	Knee	35 - 45 cm	1	6 Bags
301405	Ankle	(Shoe Size) 23 - 27 cm	1	6 Bags
301406	Shin	25 - 35 cm	1	6 Bags
301407	Thigh	45 - 55 cm	1	6 Bags

*For the elbow, measure the circumference 10 cm above the center of the elbow. For the knee, measure the circumference 12 cm above the center of the kneecap. For the ankle, use shoe size as a guide.



FACILIAID KNEE WARM SUPPORT

[Indication] ● Knee support with heat retention

3D Structure



The 3D structure is knitted to fit the leg shape. It provides equal pressure to the leg and does not obstruct blood circulation.

Active Resin Stays



The active resin stays on the both sides provide a secure fit without curling, twisting or slipping.

Heat Retention




Raised fabric is used on the inside, to trap a layer of warm air to heat the knee.

For both right and left knees

Item Number	Size	Thigh Circumference
F303101	S	40 - 43 cm
F303102	M	43 - 46 cm
F303103	L	46 - 49 cm
F303104	XL	49 - 52 cm
F303105	2XL	52 - 55 cm

● Choosing a size



Choose the size based on thigh circumference. Measure the circumference of the thigh 16 cm above the center of the kneecap in a standing position. If the measurement falls between two sizes, choose the larger one.



FACILIAID ELBOW


[Indication] ● Elbow support

- A semi-hard pad applies moderate compression on the enthesis.
- Easy to wear in one hand. Just put the arm through the ring-shaped body and fasten a hook-and-loop strap.
- The end of the strap is thicker to prevent the strap from coming out of the buckle. One size fits all.
- A washable body allows the brace to be used repeatedly and hygienically.

How to Apply



● Measuring the arm



Measure the circumference of the forearm at its widest point in a relaxed state.

For both right and left elbows

Item Number	Size	Forearm Circumference	Qty per Box
302601	One Size	19 - 31 cm	5 Bags

FACILIAID SHOULDER

FACILIAID WRIST

[Indication] ● Shoulder joint support

[Indication] ● Wrist support



Comfortable Material

Thin, well-fitting material supports the shoulder joint.



Adjustable Compression

You can adjust the support level by the attaching/detaching mechanism.

Straps on the brace and a shoulder flap can be used to adjust the amount of support.



Lifting the arm when fastening the shoulder flap, leads to a reduction in stress placed on the shoulder joint.

Heat Retention

The brace wraps around the shoulder to keep the shoulder joint warm.



A pocket warmer can be inserted between the shoulder flap and the main body.
*If a pocket warmer is used, follow the pocket warmer instructions and use it safely.




[Example]
Wearing the brace on the right shoulder



For both right and left shoulders

Item Number	Size	Shoulder Circumference
F302301	S	34 - 41 cm
F302302	M	41 - 48 cm
F302303	L	48 - 55 cm

● Choosing a size



Choose the size based on the circumference at the base of the shoulder. Measure the circumference passing over the acromion and armpit in a relaxed state. If the measurement falls between two sizes, choose the larger one.

Aluminum Stays



Palm side stay

The stay is curved to improve fit and maintain a functional position without preventing the movement of the CM/MP joints of the thumb.



Stay on the back side of the hand

Support is improved by positioning stays on both sides of the hand.



Can be worn on both right and left wrist

Both stays are removable, and the stays can be swapped over to allow the brace to be worn on either wrist.

Easy Application



Step1



Step2



The brace is put on by passing the hand through the brace and fastening the strap, so it can be put on using one hand.

The stopper on the strap prevents the strap from coming off, and makes it easier to fasten.

Breathable Material



The brace is made of soft 3D mesh fabric that fits well, breathable, and comfortable to wear.

Washable


The aluminum stays have a rust preventing finish that allows the whole brace to be washed.



For both right and left wrists

Item Number	Size	Palm Circumference
F303301	S	16 - 18 cm
F303302	M	18 - 20 cm
F303303	L	20 - 22 cm


● Choosing a size



Choose the size based on palm circumference. Measure the palm circumference around the base of the fingers, from the index finger to the little finger. If the measurement falls between two sizes, choose the larger one.

Inserting the Stays

When delivered, the stays are not inserted into the brace. Insert them in the palm side and back side of the hand according to which wrist the brace will be worn on.



Insert the stays in accordance with the enclosed instructions when giving the brace to a patient.

Check the "Palm side" and "Back side" marks printed on the stays when inserting them.

Insert the stays into the palm side and the back side of the hand.


FACILIAID THUMB

[Indication] ● Support and stabilization of the thumb CM joint

Stabilization of Thumb CM Joint



Stabilize the thumb by restricting the range of motion to reduce pain.



A built-in resin stays in the thumb part provide moderate support without interfering movement of the end of the thumb.



The movable thumb strap enables adjusting compression.

Easy Application



The brace can easily be put on using one hand by passing the hand through the brace and fastening the strap.



The plastic parts inside the wrist strap make it easy to fasten.


Thin Material

The brace has a compact shape with thin and non-bulky fabric.

Washable

The whole brace can be washed when it gets dirty.

● Choosing a size



Choose the size based on the wrist circumference at the narrowest part.
If the measurement falls between two sizes, choose the larger one.

For both right and left thumbs

Item Number	Size	Wrist Circumference	Thumb Circumference
F303401	S	12.5 - 14.5 cm	5.0 - 6.0 cm
F303402	M	14.5 - 16.5 cm	5.5 - 6.5 cm
F303403	L	16.5 - 19.0 cm	6.0 - 7.5 cm

LINE UP

2

BANDAGE PRODUCTS

- ▶ Adhesive Elastic Bandage
- ▶ Elastic Bandage
- ▶ Net & Tubular Bandage



ANKLE
KNEE
BACK
RIB
NECK
OTHERS
ADHESIVE ELASTIC BANDAGE
ELASTIC BANDAGE
NET BANDAGE/
TUBULAR BANDAGE
ALUMINIUM SPLINTS
SOFT SPLINTS
OTHERS

ANKLE
KNEE
BACK
RIB
NECK
OTHERS
ADHESIVE ELASTIC BANDAGE
ELASTIC BANDAGE
NET BANDAGE/
TUBULAR BANDAGE
ALUMINIUM SPLINTS
SOFT SPLINTS
OTHERS




Moderate Elasticity
BENEFACT Adhesive Elastic Bandage

[Primary uses]

- Stabilization after sprain or dislocation
- Securing tubes such as catheters or drains
- Postoperative compression and stabilization

- The bandage has moderate elasticity and adhesion to fit active joints.
- Acrylic adhesive that is gentle to the skin is used.
- Two colors are available: White, Beige.

Item Number	Type		Specifications (Width × Actual length)	Qty per Box
	White	Beige		
4101	4111	1	1.25 cm × 5 m	24 Rolls
4102	4112	2	2.50 cm × 5 m	12 Rolls
4103	4113	3	5.00 cm × 5 m	6 Rolls
4104	4114	4	7.50 cm × 5 m	4 Rolls
4105	4115	5	10.00 cm × 5 m	3 Rolls
4106	4116	6	3.75 cm × 5 m	8 Rolls



1 cm squares on the release paper act as a cutting guide.




Soft, Non-Woven Fabric
BENELUKE Adhesive Gauze Bandage

[Primary uses]

- Securing gauze, absorbent cotton, and poultices
- Securing tubes such as catheters or drains

- Acrylic adhesive that is gentle to the skin is used.
- The adhesive is applied in stripes to improve breathability.

Item Number	Type	Specifications (Width × Actual length)	Qty per Box
17902	2	2.5 cm × 10 m	12 Rolls
17904	4	5.0 cm × 10 m	6 Rolls
17905	5	7.5 cm × 10 m	4 Rolls
17906	6	10.0 cm × 10 m	3 Rolls
17907	7	15.0 cm × 10 m	2 Rolls



▲ 1 cm squares on the release paper act as a cutting guide.




Moderate Compression and Stabilization
BENESTORK Adhesive Elastic Bandage

[Primary uses]

- Securing tubes such as catheters or drains
- Compression and stabilization after surgery, to stop bleeding, for puncture and rib fractures
- Restriction and stabilization of range of joint motion
- Taping for exercise

- Hard-twist yarn is used as a thick base cloth.
- Suitable for compression and stabilization where high elasticity is required.
- Acrylic adhesive that is gentle to the skin is used.
- The adhesive is applied in stripes to improve breathability.

Item Number	Type	Specifications (Width × Actual length)	Qty per Box
4302	2	2.50 cm × 5 m	12 Rolls
4303	3	3.75 cm × 5 m	8 Rolls
4304	4	5.00 cm × 5 m	6 Rolls
4305	5	7.50 cm × 5 m	4 Rolls
4306	6	10.00 cm × 5 m	3 Rolls



▲ 1 cm squares on the release paper act as a cutting guide.



100% Cotton Bandage
L-WEAVE Elastic Bandage

[Primary uses]

- Securing splints
- Postoperative compression and immobilization

- The extension rate is approximately 200%.
- The bandage is made from thick material, making it ideal for compression and securing splints.
- 100% cotton fabric with excellent moisture absorption.
- Sterilization by autoclave is possible.
- No fluorescent coating is used.

Item Number	Type	Specifications (Width × Extension)	Qty per Box
022	2	5.0 cm × 4.5 m	6 Rolls
023	3	7.5 cm × 4.5 m	6 Rolls
024	4	10.0 cm × 4.5 m	6 Rolls
027	7	17.5 cm × 4.5 m	6 Rolls



Ideal Elasticity and Compressing
ELTECT Elastic Bandage

[Primary use]

- Compression and stabilization for burns and after surgery

- The extension rate is approximately 250%.
- It can be used on various body parts with moderate compression and stabilization.
- Suitable for compression and stabilization for burns or after surgery.
- Highly stretchy fibers are used to provide high elasticity.

Item Number	Type	Specifications (Width × Extension)	Qty per Box
052	2	5.0 cm × 9 m	10 Rolls
053	3	7.5 cm × 9 m	10 Rolls
054	4	10.0 cm × 9 m	10 Rolls
056	6	15.0 cm × 9 m	4 Rolls



Excellent Elasticity
L-WHITE Elastic Bandage

[Primary uses]

- Securing gauze and poultices
- Protection of wound surfaces

- The extension rate is approximately 300%.
- Suitable as a dressing for everyday wounds.
- Highly stretchy fibers are used to allow the bandage to fit active joints well.
- The durable knitted fabric prevents from fraying.

Item Number	Type	Specifications (Width × Extension)	Qty per Box
012	2	5.0 cm × 9 m	10 Rolls
013	3	7.5 cm × 9 m	10 Rolls
014	4	9.0 cm × 9 m	10 Rolls



Anti-Slipping Bandage
L-WHITE-S Elastic Bandage

[Primary uses]

- Securing gauze and poultices
- Protection of wound surfaces

- The extension rate is approximately 300%.
- Suitable as a dressing for everyday wounds.
- Highly stretchy fibers are used to allow the bandage to fit active joints well.
- The woven texture provides an uneven surface that prevents from slipping.
- Comfortable, breathable fabric with high moisture absorption.

Item number	Type	Specifications (Width × Extension)	Qty per Box
032	2	5.0 cm × 9 m	10 Rolls
033	3	7.5 cm × 9 m	10 Rolls
034	4	9.0 cm × 9 m	10 Rolls



RETELATAY

Elastic Net Bandage

[Primary uses]

- Securing gauze and poultices

- Stretch net tubular bandage.
- Can be cut to any shape based on the affected area.
- The bandage is made using a unique knitting technique that prevents fraying on the cut edges.
- A variety of sizes are available, from finger size to body trunk size. Size 1 comes with an applicator useful for finger treatment.

Contains Latex Rubber



◀ The bandage can be applied without touching the affected area.

Item Number	Type	Body Part	Specification (Width x Extension)
00061	1	Finger	0.7 cm x 25 m
00062	2	Limb	2.0 cm x 25 m
00063	3	Arm, Lower Leg, Foot	2.4 cm x 25 m
00064	4	Shoulder, Thigh, Knee	3.2 cm x 25 m
00065	5	Head, Shoulder, Thigh	5.5 cm x 25 m
00066	6	Head, Rib, Waist	7.0 cm x 25 m
00067	7	Rib, Waist	8.0 cm x 25 m

ROLL-UP STOCKINET

Tubular Bandage

[Primary uses]

- Keeping limbs warm during surgery

- The bandages come pre-rolled into a donut shape. They can be placed over hands or legs and rolled down for quick and accurate application.
- Made from high-quality pure cotton with high moisture absorption and a soft texture.
- Six different types of bandages allow all four limbs to be covered. The different types are color coded for easy differentiation.
- Sterilization by autoclave or EOG is possible.



Edges are color coded according to length ▶

- 150 cm
- 120 cm
- 90 cm



Item Number	Type	Coverage	Specifications (Width x Length)	Qty per Box	Edge Color
110030	3	Upper Limbs	7.6 x 90 cm	20	White
110041	4A		10.2 x 90 cm	20	White
110042	4B	Lower Limbs	10.2 x 120 cm	20	Green
110043	4C		10.2 x 150 cm	20	Blue
110052	5B		13.7 x 120 cm	20	Green
110053	5C		13.7 x 150 cm	20	Blue


CASTING PRODUCTS

- ▶ Aluminum Splints
- ▶ Soft Splints
- ▶ Others



ANKLE
KNEE
BACK
RIB
NECK
OTHERS
ADHESIVE ELASTIC BANDAGE
ELASTIC BANDAGE
NET BANDAGE/TUBULAR BANDAGE
ALUMINIUM SPLINTS
SOFT SPLINTS
OTHERS

BRACES & SUPPORTS
Ankle
Knee
Back
Rib
Neck
Others
NURSING CARE PRODUCTS
Adhesive Elastic Bandage
Elastic Bandage
Net Bandage/
Tubular Bandage
CASTING PRODUCTS
Alumi Splints
Soft Splints
Others



ALUMINUM SPLINTS

[Indication]

- Splinting of fractures and contusions

- Efficient splints made of a combination of aluminum alloy plates and urethane foam cushioning.
- Impact to the affected area is absorbed by the comfortable cushioning material.
- Suitable for small body parts such as fingers or arms. Available in 12 different shapes and sizes.
- Can be bent to fit the affected area. Provides secure protection and stabilization.



ALUMINUM SPLINTS

Item Number	Type	Specifications (Depth x Width x Length)	Qty per Box
5001	8A (Frog Shape)	1.0 x 57 mm	6
5002	8B (Frog Shape)	1.0 x 70 mm	6
5005	10A (Plate Type)	1.5 x 13 x 197 mm	12
5006	10B (Plate Type)	1.5 x 13 x 400 mm	6
5009	12A (Plate Type)	1.5 x 18 x 197 mm	12
5010	12B (Plate Type)	1.5 x 18 x 400 mm	6
5012	13B (Plate Type)	1.5 x 25 x 400 mm	6
5014	2B (Plate Type)	1.5 x 50 x 400 mm	6
5015	15 (Spoon Shape)	1.0 x 20 x 40 mm	12
5016	16 (Spoon Shape)	1.0 x 20 x 67 mm	12
5017	17 (Spoon Shape)	1.0 x 20 x 95 mm	6
5018	18 (Spoon Shape)	1.0 x 20 x 150 mm	6



SOFT SPLINTS

[Indication]

- Splinting for fractures and sprains

- User-friendly design.
- Can be bent where necessary.
- Made from galvanized iron wire covered with cushioning for a comfortable fit. No cover is required.
- Comfortable material with good moisture absorption is used where the splint touches the skin.
- Available in 8 different sizes to fit the affected area.

Item Number	Type	Specifications (Width x Length)	Qty per Box
5499	5S	4.5 x 21.5 cm	6
5500	4S	4.5 x 32.0 cm	6
5501	3S	6.0 x 41.0 cm	6
5502	SS	7.0 x 51.0 cm	6
5503	S	8.0 x 63.0 cm	6
5504	M	9.0 x 72.0 cm	6
5505	L	10.0 x 83.0 cm	6
5506	LL	12.0 x 113.0 cm	6



ARM SLING

[Primary use]

- Holding and protection of upper limb

- An adjustable pad to reduce pressures from your shoulder for enhanced comfort.
- Made with a thin, breathable material for extra comfort.
- Designed to fit the elbow shape when wearing the sling.
- The thin, reversible buckle reduces bulk. The sling can be worn on either arm.
- The elbow and wrist areas are reinforced with a dual structure and straps to reduce strain.

For both right and left arms

Item Number	Size	Forearm Length
F601901	L	20 - 30 cm
F601902	M	- 20 cm

*Choose the size based on forearm length. Measure the forearm length from the olecranon to the wrist. If the measurement falls between two sizes, choose the larger one. Medium size is suitable for child.



CAST BOOT

[Primary use]

- Preventing casts or splints for ankle joints from becoming wet or stained

- The ship bottom shaped sole is easy to walk in and stable. Striped finishing prevents slipping.
- Lightweight, triple layered structure with excellent cushioning.
- The hook-and-loop fastener makes the boot easy to put on and take off.
- The waterproof body prevents the cast from getting wet.

For both right and left feet

Item Number	Size	Shoe Size	Specifications (Width x Depth)
6014	S	19.0 - 22.0 cm	11.0 - 24.5 cm
6015	M	22.0 - 25.0 cm	12.5 - 26.5 cm
6016	L	25.0 - 28.0 cm	13.5 - 31.0 cm

*We recommend trying the shoes on when choosing a size.